

Hel shaqo! Eeg talaabada ugu horeysa ee caafimaad oo sanbab oo dhaama

1 Isku xidhxidh spirometer oo geli baytariyada

2 Soo dejiso ablikeeshinka oo samayso xisaab



Ablikeeshinka waxaa laga heleyaa iOS iyo Android. Iskaan garee QR-ka koodhka si aad ablikeeshinka u hesho

3 Aqbal casuumada xarunta caafimaadka

Marka aad aqbalayso casuumada xarunta caafimaadka, waxa aad ogolaanaysaa inaad siiso xogta shaybaadhkaaga xirfadlaha daryeelka caafimaadka xarunta caafimaadka

4 Hel qalabka AsthmaTuner

Miyay kugu bateen xasuusintu? Qalabka xasuusinta ee AsthmaTuner ayaa kuu xasuusin doonaa si xiiso leh inaad qabato shaqooyinka oo aad ku dhamaystirto ablikeeshinka

Widget AsthmaTuner



Qalabyadaan waxay tusaan macluumaadka ablikeeshinka iyada oo aanad furin ablikeeshinka laftiisa. Waa iska caadi in la yeesho qalabka cimilada, waxa uu si toos ah kuugu tusayaa saadaasha cimilada shaashada taleefankaaga.



Qalabyadeenu waxay tusaan sida geedku yahay, waxa uu kugu xasuusiyaan sida ugu wanaagsan in shaybaadhkaaga aad marto ama aad samayso wax inta aad gasho ablikeeshinka AsthmaTuner.